IN-BIOS[9,5]000 2022

#### Basic file formats

Arvind Sundaram Oct 18, 2022

Norwegian Sequencing Centre OUS, Ullevål, Oslo

#### File formats

- \* FASTA sequence data
- \* FASTQ sequence data including quality (SRA)
- \* GTF/GFF, GFF3 Gene feature & transfer format
- \* SAM/BAM aligned/mapped data (CRAM)
- \* VCF/BCF variant calling

### FASTA

- \* Represents a sequence DNA/RNA/protein
- \* Two lines:
  - \* header starts '>'
  - Second line contains the sequence

header

>NM\_001199335.1 Danio rerio toll-like receptor 21 (tlr21), mRNA CTAGACGAAGCTCTCACAGAGTGAAAGGCTTTGTAATGTTGGACTCAGTTTATAGGAATT ACATTTAAAAGTGGAACTGAAAAACATGGCACACTCTGCATGCCACAAACTGATACTAAA GGCCACATTCATCTGTCTCATAAAACTTGCCTGCAGCTACAGTTTCAGAAGTTGCATAGA GATCCCAGATTCTAATCATACAATCTTTACATGTGTTAAAAGTTATGAACGAGACATAAC TGCGATTGTGAGTGATGTATTTCCCACTGCATTAAATCTTACAATCTCTCACA

Genomes Transcriptomes

### FASTQ

- \* Represents a sequence along with quality information
- \* Four lines:
  - header starts '@'
  - \* Second line contains the sequence
  - \* Quality header starts with '+'
  - \* Fourth line contains the quality for each base in the second line




# FASTQ quality encoding

```
!"#$%&'()*+,-./0123456789:;<=>?@ABCDEFGHIJKLMNOPQRSTUVWXYZ[\]^_`abcdefghijklmnopqrstuvwxyz{|}~
33
                                         126
                               104
Phred+33, raw reads typically (0, 40)
S - Sanger
     Solexa+64, raw reads typically (-5, 40)
X - Solexa
I - Illumina 1.3+ Phred+64, raw reads typically (0, 40)
J - Illumina 1.5+ Phred+64, raw reads typically (3, 41)
 with 0=unused, 1=unused, 2=Read Segment Quality Control Indicator (bold)
 (Note: See discussion above).
L - Illumina 1.8+ Phred+33, raw reads typically (0, 41)
```

### GFF/GTF


- \* GFF: General Feature Format
  - Current format is based on version 3 and called GFF3
  - format consists of one line per feature
  - \* each containing 9 columns of data, plus optional track definition lines.
- \* GTF: General Transfer Format is identical to GFF version 2.5.

#### GFF3



#### BED format

- Browser Extensible Data
- \* 3-12 fields; but generally contains 6 or more



There are far more annotation formats but will stop here!!!

## Chromosome coordinate systems

\* 0-based

**ACTGACTG** 

01234567

\* 1-based

**ACTGACTG** 

12345678

To represent the TGAC:

0-based inclusive: 2-5

1-based inclusive: 3-6

1-based exclusive: 3-7

Ensembl: 1-based inclusive

UCSC: 0-based start

1-based end

1-based for display

SAM, VCF, GFF, GTF: 1-based

BAM, BCF: 0-based

BED: 0-based exclusive

Most tools are aware of these differences

### SAM/BAM

- \* SAM Sequence Alignment/Map format
- \* BAM same as/similar to SAM but in binary format
  - Most softwares and tools prefer BAM

\* Header

(access using samtools view -H)

Alignment records

(access using samtools view; use -h to include Header)

### SAM/BAM: header

```
SO:coordinate ← sort order
@HD
    VN:1.0GO:none
@SQ
     SN:chrM LN:16571
@SO
     SN:chr1 LN:249250621
                                       Reference sequence names with length information
@SO
     SN:chrX
             LN:155270560
@SQ
     SN:chrY
             LN:59373566
                                                       Read group with platform,
    ID:86-191PL:ILLUMINA LB:IL500 SM:86-191-1
@RG
                                                      library and sample information
@RG
    ID:BsK010PL:ILLUMINA LB:IL501 SM:BsK010-1
                                                         Program analysis history
@RG
    ID:SDH023PL:ILLUMINA LB:IL508 SM:SDH023
@PG
     ID:GATK IndelRealigner VN:2.0-39-gd091f72CL:knownAlleles=[]
targetIntervals=tmp.intervals.list LODThresholdForCleaning=5.0
consensusDeterminationModel=USE READS entropyThreshold=0.15
maxReadsInMemory=150000 maxIsizeForMovement=3000 maxPositionalMoveAllowed=200
maxConsensuses=30 maxReadsForConsensuses=120 maxReadsForRealignment=20000
noOriginalAlignmentTags=false nWayOut=null generate_nWayOut_md5s=false
check early=false noPGTag=false keepPGTags=false indelsFileForDebugging=null
statisticsFileForDebugging=null SNPsFileForDebugging=null
   ID:bwaPN:bwaVN:0.6.2-r126
aPG
```

http://samtools.github.io/hts-specs/SAMv1.pdf

# SAM/BAM: alignment records

1
HW-ST605:127:B0568ABXX:2:1201:10933:3739
2 3 4 5 6 7 8 9
HW-ST605:127:B0568ABXX:2:1201:10933:3739
147 chr1 27675 60 101M = 27588 -188

- 10 TCATTTTATGGCCCCTTCTTCCTATATCTGGTAGCTTTTAAATGATGACCATGTAGATAATCTTTATTGTCCCTCTTTCAGCAGAC
- 11 =7;:;<=??<=BCCEFFEJFCEGGEFFDF?E@E@ADCACB>CCDCBACDCDDDAB@@BCADDCBC@BCBB8@ABCCCDCBDA@>:/

Col	Field	Type	Regexp/Range	<b>Brief Description</b>
1	QNAM	String	[!-?A-~]{1,254}	Query template NAME
2	FLAG	Int	$[0,2^{16}-1]$	bitwise FLAG
3	RNAME	String	\* [!-()+-<>-~][!-~]*	Reference sequence NAME
4	POS	Int	$[0,2^{31}-1]$	1-based leftmost mapping POSition
5	MAPQ	Int	[0,28-1]	MAPping Quality
6	CIGAR	String	\* ([0 <b>-</b> 9]+	CIGAR string
7	RNEXT	String	\*   =   [!-()+-<>-~][!-~]*	Ref. name of the mate/next read
8	PNEXT	Int	$[0,2^{31}-1]$	Position of the mate/next read
9	TLEN	Int	$[-2^{31}+1,2^{31}-1]$	observed Template LENght
10	SEQ	String	$\ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ $	sequence SEQuence
11	QUAL	String	[!-~]+	ACSII of Phred-scaled base QUALity+33

http://samtools.github.io/hts-specs/SAMv1.pdf

### SAM/BAM

\* To check FLAG info: <a href="https://broadinstitute.github.io/picard/explain-flags.html">https://broadinstitute.github.io/picard/explain-flags.html</a>

\* BAM/SAM files are produced by alignment/mapping softwares and other tools. You will hear more about these in Tim Hughes's *Algorithms* lecture

- QC of BAM/SAM files Picard, Qualimap
- Manipulating and extracting information from varied file formats -SAMtools, BEDtools, BCFtools

### VCF/BCF

- \* Variant call format
- \* BCF same as/similar to VCF but in binary format
  - \* Most softwares and tools prefer BCF or VCF in zip format
- header and records

```
##fileformat=VCFv4.2
##source=myImputationProgramV3.1
##reference=file:///seq/references/1000GenomesPilot-NCBI36.fasta
##contig=<ID=20,length=62435964,assembly=B36,md5=xxxx,species="Homo sapiens",taxonomy=x>
...
##INFO=<ID=NS,Number=1,Type=Integer,Description="Number of Samples With Data">
##INFO=<ID=DP,Number=1,Type=Integer,Description="Total Depth">
...
##FILTER=<ID=q10,Description="Quality below 10">
...
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
...
##CHROM POS ID REF ALT QUAL FILTER INFO FORMAT NA00001 NA00002 NA00003
```

https://samtools.github.io/hts-specs/VCFv4.2.pdf

### VCF/BCF: records

#CHROM POS ID REF ALT QUAL FILTER INFO FORMAT NA00001 20 14370 rs6054257 G A 29 PASS NS=3;DP=14;AF=0.5;DB;H2 GT:GQ:DP:HQ 0 0:48:1:51,51

20 17330 . T A 3 q10 NS=3;DP=11;AF=0.017 GT:GQ:DP:HQ 0 0:49:3:58,50

20 1110696 rs6040355 A G,T 67 PASS NS=2;DP=10;AF=0.333,0.667;AA=T;DB GT:GQ:DP:HQ 1 2:21:6:23,27

	#CHROM	POS	ID	REF	ALT	QUAL	FILTER	INFO	FORMAT	SAMPLES
SNP	20	3		С	G		PASS	DP=100		
Deletion	20	2		TC	С		PASS	DP=100		
Insertion	20	2		TC	TCA		PASS	DP=100		
Alleles	20	2		TC	TG,T		PASS	DP=100		
Alleles	20	2		TCG	TG,T,TCAG	•	PASS	DP=100		

#### + Structural variants

http://vcftools.sourceforge.net/VCF-poster.pdf

# Qualimap

- \* GUI, command line based
  - examines sequencing alignment data in SAM/ BAM format
- \* requires Java, R
- \* Uses:
  - FastQC, Picard, SAMTools, NOISeq, Repitools and JProfiler

- \* BAM QC
- \* RNA-seq QC
- Counts QC
- \* Multi-sample BAM QC

### SAMtools/BCFtools

- \* SAMtools
  - manipulating SAM/BAM
- \* BCFtools
  - manipulating VCF/BCF
- \* Part of HTSlib tool kit

```
-- indexing
            index/extract FASTA
   faidx
            index alignment
   index
-- editing
             recalculate MD/NM tags and '=' bases
   calmd
   fixmate
             fix mate information
   reheader replace BAM header
             remove PCR duplicates
   rmdup
   targetcut cut fosmid regions (for fosmid pool only)
-- file operations
   bamshuf shuffle and group alignments by name
           concatenate BAMs
   cat
             merge sorted alignments
   merge
   mpileup multi-way pileup
           sort alignment file
   sort
   split
           splits a file by read group
             converts a BAM to a FASTQ
   bam2fq
-- stats
   bedcov
             read depth per BED region
            compute the depth
   depth
   flagstat simple stats
   idxstats BAM index stats
            phase heterozygotes
   phase
           generate stats (former bamcheck)
   stats
-- viewing
   flags
           explain BAM flags
            text alignment viewer
   tview
            SAM<->BAM<->CRAM conversion
   view
```

#### **BED**tools

- \* Awesome tool to compare dataset in BED/GFF/VCF and BAM formats
- Very well documented and Highly recommended!!

\* Example:

http://bedtools.readthedocs.io/en/latest/content/tools/intersect.html



# IGV - Integrative Genomics Viewer


\* View BAM, VCF, GTF files and more

Ensembl browser

http://ensembl.org/

\* UCSC browser

http://genome.ucsc.edu/



http://software.broadinstitute.org/software/igv/

# Further reading

- \* GFF3, GTF <a href="http://gmod.org/wiki/GFF3">http://gmod.org/wiki/GFF3</a>
- \* SAM/BAM format <a href="http://samtools.github.io/hts-specs/SAMv1.pdf">http://samtools.github.io/hts-specs/SAMv1.pdf</a>
- VCF/BCF format <a href="https://samtools.github.io/hts-specs/VCFv4.2.pdf">https://samtools.github.io/hts-specs/VCFv4.2.pdf</a>
  - <a href="http://vcftools.sourceforge.net/VCF-poster.pdf">http://vcftools.sourceforge.net/VCF-poster.pdf</a>
- \* IGV <a href="http://software.broadinstitute.org/software/igv/">http://software.broadinstitute.org/software/igv/</a>

#### Databases

- Ensembl <a href="http://ensembl.org/">http://ensembl.org/</a>
- \* UCSC <a href="http://genome.ucsc.edu/">http://genome.ucsc.edu/</a>
- \* NCBI <a href="http://ncbi.nlm.nih.gov/">http://ncbi.nlm.nih.gov/</a>
  - \* RefSeq https://www.ncbi.nlm.nih.gov/refseq/

There are too many databases containing specific information!!!